


Coherent X-ray Imaging Advanced Light Source, BL 9.0.1



Introduction

Diffractive Imaging (Yeast Cells, 3D, Aerogels, Ultrafast)

Massively Parallel X-ray Holography

Femtosecond Time-Delay Holography


Serial Crystallography and Particle Injection



Motivation

- Obtain high resolution structural information about biological and material science specimens
- Develop experimental techniques optimized for upcoming coherent X-ray lasers

Beamline 9.0.1




- A pinhole selects a monochromatic coherent x-ray beam
- A secondary pinhole cleans the beam
- X-rays: 520-750 eV (1500 eV upgrade in progress)
- Sample can be rotated perpendicular to x-rays
- Area detector collects the diffraction pattern *intensity*



Diffraction Imaging

Diffraction imaging is a form of microscopy whereby a lens is replaced by an array detector and a computer. The computer recovers the phases missing in the recorded intensity patterns.




Proposed by D. Sayre, first demonstrated by J. Miao et al, in 1999



Overview Slide: Diffractive Imaging


Experimental

We can collect diffraction patterns at many orientations at the Advanced Light Source (wavelength: 1.6 nm)



We solve the phase problem
For 10^8 reflections



**Stony Brook developed
endstation for cryo imaging**



Yeast cell imaged at 40 nm resolution, D.Shapiro et al., *Proc. Nat. Acc. Sci. U.S.A.* 102, 15343-15346, (2005).




Test object imaged in 3D



Chapman et al. J. Opt. Soc. Am A 23, 1179 (2006)



2D view





Algorithms

Shrinkwrap algorithm, S. Marchesini et al, *Phys Rev B* 68, 140101(R) (2003)



Ta₂O₅ Nanofoam



Structural properties of Aerogels revealed. Barty et al. submitted





Diffractive Imaging of Yeast cells

Motivation: obtain high resolution 3D images of cells


Yeast cell: 2.5 micron thick,
unstained freeze-dried, at 750 eV
Total dose $\sim 10^8$ Gray (room temperature)
Oversampling is about 5 in each dimension



Gatan 630 cryo holder keeps cells frozen



David Shapiro, Stony Brook, now at ALS





Consistency of reconstructions
quantifies the resolution: 30 nm.

3D Diffraction-imaging


- Complete coverage of reciprocal space by sample rotation
- Use a true 3D object that can be well-characterized by independent means

Diffraction patterns collected at multiple orientations



FFT

Autocorrelation obtained at multiple orientations




$\pm 70^\circ$ range
 1° angular spacing


Multiple views
(movie)

3D Reconstruction Is Achieved by Fourier Synthesis


One diffraction pattern gives information on the Ewald sphere in reciprocal space




Rotating a sample about one axis allows to fill the data in Fourier space



Full 3D Image Reconstructions From Experimental X-ray Data



Voxel Data cube 1024^3 elements



A. Barty



Reconstruction



Chapman et al. JOSAA 23, 1179 (2006)

Size	Memory needed	
	Single Precision	Double Precision
2563	336 MB	592 MB
5123	2.6 GB	4.7 Gb
10243	22 GB	38 GB
20483	176 GB	304 GB

Size	32-CPU G5 cluster speed	
	3D FFT timing	Processing time*
2563	73 msec	10 mins
5123	850 msec	1.5 hrs
10243	7.9 sec	14 hrs
20483	??	??


*2000 iterations, 2FFTs per iteration plus other floating point operations code: C with mvapich, dist_fft (Altivec-optimized MPI FFT) see http://images.apple.com/acg/pdf/20040827_GigaFFT.pdf

2D Single-view Images Have Depth Information





Point spread function (3D)

Propagation of wavefield




True 2D Projection Images Is Formed From a Central Section of the 3D Diffraction Data



Infinite depth of focus

A true projection image is obtained from a **plane central section** of the 3D diffraction data. Data must be collected at many object orientations to achieve this




Coherent X-ray Diffraction Imaging of Aerogel Structures Compared With USAXS and TEM

Method	Rods radius (nm)	Rods length (nm)	Isotropic?
Diffraction Image	<13.7	44	no
USAXS	6.8	40	no ?
TEM	3-6	10-40	Can't be determined

Motivation



Aerogel strength/weight ratio depends on structural properties. We find a structure formed of blobs and rods, indicating pathways to make stronger and lighter materials



M. Wang (LLNL)

TEM results based on survey of 3-4 cells;
Diffraction based on ~300 cells.

USAXS





A. Barty (LLNL), submitted


Ultrafast Coherent Diffractive Imaging

Motivation: To study ultrafast dynamics and overcome radiation damage
We applied the techniques developed at ALS at the first X-FEL source

ALS, $\lambda=1.6$ nm





SEM



X-rays


Chapman et al. JOSAA 23, 1179 (2006)

FLASH, $\lambda=32$ nm



SEM



1 micron




Single pulse
FEL

Chapman Et Al. Nature Physics 2, 839 - 843 (2006)

Image Reconstructed From an Ultrafast FEL Diffraction Pattern




NA = 0.26. This is the highest numerical aperture coherent X-ray diffraction pattern ever reconstructed. Pixel size is 60 nm.



Reconstructed Image – achieved diffraction limited resolution!
Wavelength = 32 nm



Overview Slide: Massively Parallel Holography With Coded Apertures





We improved the efficiency of holography by a factor of 2500

Signal to noise improvement

Fourier Transform Holography

The Fourier transform of the diffraction pattern provides the autocorrelation of the object, and a reference point generates a hologram in the autocorrelation map




- The smaller the dot, the higher the resolution
- But the smaller the dot, the dimmer the hologram

Coded Apertures Overcome Resolution Vs Brightness Limitations


One point creates a hologram, many points create overlapping holograms:
Like a pinhole camera with many pinholes.

The “magic trick”:
An extended object with
point-like autocorrelation
(uniformly redundant array)




Coded apertures are used in many applications
Where lenses do not exist to improve signal

Pinhole camera




γ -ray astronomy




Medical imaging



Homeland security





spectroscopy





Ultrafast, Ultrabright, Ultra High Resolution X-ray Holography With Coded Apertures



$\lambda=2.2$ nm Resolution=43 nm, SNR X ~20



$\Delta t=10$ fs, $\lambda=13.5$ nm, Res~100 nm




Overview Slide: Femtosecond Time Delay Holography

Motivation: to study the ultrafast dynamics of particles under intense x-ray fields



Scattering from a dusty mirror corroborates Livermore code for plasma dynamics




Scattering of 97 nm spheres ALS



Size selection measured at ALS




Time delay holography, FLASH



Plasma dynamics measured at FLASH

Density (g/cm^3)




Successful experiments at the ALS have validated sample purification strategy and provided insight into data FLASH FEL expts.

Mechanism of X-ray Induced Damage


Motivation: To study the ultrafast dynamics of particles under intense x-ray fields

Ionization leads to charge redistribution and explosion

- Some high-energy electron escape the molecule and charge it positively
- Remaining free electrons equilibrate and redistribute




Trapped electrons shield the core, forming a 2-layer configuration



Selecting Identical Particles


Gas phase mobility separation for preparing size-monodisperse polystyrene nanoparticles



Polystyrene Sphere Scattering Measured at the Advanced Light Source Validated <5% Size Variation


Successful experiments at the ALS have validated sample purification strategy and provided insight into data expected at upcoming VUV-FEL expts

Scattering of 97 nm spheres




Wavelength = 1.65 nm
See 7th minima distinctly

Line-out Average of Diffraction Patterns




Easily differentiate particles with $D = 97$ and 102 nm, model shows <3.5% variation


Particle Explosion Experiments Were Performed at FLASH on Latex Particles on Membranes




Mounted on piezo x-y stage to move each window into beam





Latex particles on membrane




LLNL VUV Hydrodynamic Code Shows That Latex Spheres Start Exploding in ~ 2 Ps



Scattering From Balls Demonstrates That They Retain Their Shape Throughout the Duration of the Pulse





$\lambda = 32 \text{ nm}$



$\lambda = 13 \text{ nm}$

q (μm^{-1})

We Applied Femtosecond Time-delay Holography to Study Explosion Process





The hologram is reflected onto an area detector


— Prompt diffraction
— Delayed diffraction

Time delay $2l/c$

First Demonstration of Time-delay Holography With 3 Fs Time Resolution Indicates the Particle Explosion




The “dusty mirror” experiment




Single shot ultrafast time-delay X-ray hologram, with 300 fs delay

The Explosion Is in Good Agreement With Our Hydrodynamic Model

Experiments





Calculation





- The structure factor narrows, showing the particle exploding
- The lower resolution shape of the explosion is different than expected
- This is the first high-resolution observation of particle explosions

We Interferometrically Measure the Change in Optical Density of the Particle at Short Delays




Injection Strategies for Imaging in Vacuo


Charge-reduced electrospray of nanoparticles




Calculated particle trajectories



Injection system successfully used with FELs




Mike Bogan, et al. Nanoletters, in press



Validated
Size selection


Triggered Aerojet injection tested at ALS



Weierstall, et al. NIM A in press

Serial Crystallography: Laser Alignment of Proteins


Motivation: To Image Proteins That Can't Be Crystallized




- Protein-in-water droplets injected into vacuum and rapidly cool by evaporation to form vitreous ice
 - A polarized laser beam aligns the protein due to induced polarizability of the molecule
 - X-ray diffraction is accumulated for many thousands of molecules
 - When enough signal is accrued, polarization is rotated and diffraction recorded at a different angle
 - 3D molecule image is obtained by phase retrieval (lensless imaging)
- If alignment is good, one could use a standard synchrotron:
- Aligned molecules briefly pass through X-ray beam (exposure less than damage threshold)

Diffraction rings

Worlds smallest
protein nanocrystals
(Photosystem 1)





Weierstall et al. NIM A
in press.




Shapiro et al summer 07
ALS 9.0.1

Laser alignment of injected proteins

To solve proteins which can't be crystallized.
Method: Spray proteins, align by laser and image.



Computational Lenses




ALS organized the first and last of a series of international conferences on coherent x-ray imaging



We are leading the development of reconstruction methods



We Improved the Speed and Reliability of Reconstruction Methods



Complex Lena

2D Step size optimization
4D Step size optimization

Collaborations

LBL: Malcolm Howells, Janos Kirz, David Shapiro, Stefano Marchesini

Stony Brook: Xiaojing Huang, Chris Jacobsen, Janos Kirz, Enju Lima Huijie
Miao Aaron M. Neiman, Johanna Nelson, David Sayre,
Jan Steinbrener, Andrew Stewart

SLAC: Sebastien Boutet

LLNL: Anton Barty, Henry Chapman, Stefan Hau-Riege,

ASU: John Spence, Uwe Weierstall, Bruce Doak, Dan Deponete

