



An Approach to High Resolution *Ex-Situ* NMR


NMR spectrum of *trans*-2-pentenal under ordinary high resolution conditions. Peaks corresponding to the five non-equivalent hydrogen nuclei are observed.




Spectrum with sample placed in non-uniform field outside the rf coil; all chemical information is lost.




Specially designed rf pulse train (inset, rectangles) "refocuses" nuclei in nonuniform magnetic field; high resolution spectrum is recovered, from emitted signal (horns).




In traditional NMR, the sample (blue) is placed in the uniform magnetic field inside the RF coil. **Ex-situ** conditions were simulated by placing the sample outside of this region.



Potential Future Applications



Scanning micro-tip
NMR analysis



Ex-situ NMR
"in the field"



Ex-situ MRI
"in the doctor's office"